



### Ikegami MKC-500HD Full-HD Medical Camera




### Ikegami Medical Monitor



### Roland VC-HD50 converter



### Roland FI HDD recorder



necessary to convert from M2T to MOV by ClipWrap2(\$49.99)  
<http://www.divergentmedia.com/clipwrap>

### Apple Final Cut Pro / iMovie



### DVD authoring by DVD Studio Pro or iDVD


## Ikegami Medical Monitor



## Ikegami MKC-500HD Full-HD Medical Camera





## Ikegami GFS-V10 Flash Memory Recorder



necessary to convert from MXF to MOV by Calibrated MXF{Q}Import (\$249.95)  
<http://www.calibratedsoftware.com>

## Apple Final Cut Pro



DVD authoring by DVD Studio Pro

Capacity	Recording Time				
	MPEG-SD30 (720x480)	MPEG-SD40 (720x480)	MPEG-SD50 (720x480)	MPEG-HD50 (1920x1080)	MPEG-HD100 (1920x1080)
16GB	43 min.	40 min.	36 min.	36 min.	18 min.
32GB	86 min.	80 min.	72 min.	72 min.	36 min.
64GB	172 min.	160 min.	144 min.	144 min.	72 min.
128GB	344 min.	320 min.	288 min.	288 min.	144 min.

GFAK media



CF Adapter